

Performance comparison of batch and multicolumn protein A capture processes


Summary

- Batch, 2C-PCC (CaptureSMB), 3C-PCC and 4C-PCC Protein A capture processes were compared in a numerical optimization study using an experimentally validated shrinking core model for mAb capture
- Performance parameters for comparison: Productivity (Throughput) and Capacity Utilization (Resin costs)
- Outcome 1: Multicolumn processes ensure high capacity utilization, low buffer consumption
- Outcome 2: 2C-PCC outperforms 3C-PCC and 4C-PCC
- Additional advantage of 2C-PCC: Minimal multi-column hardware configuration

Process Description


- Single column Batch Capture: A large fraction of the affinity resin remains unused
- Used capacity typically 40-50% of static binding capacity (SBC).
- Sequential loading processes: By splitting the bed into two columns, the capacity utilization is maximized.
- Used capacity typically > 90% of SBC
- Faster loading flow rates can be used

Multicolumn processes complexity



Chromatography model


- Shrinking core model:



- Resin: MabSelect Sure
- Harvest titer 2.5 g/L // 5.0 g/L
- 10 cm bed height per column
- CIP time 30 min
- mAb: IgG1

Model calibration

- The model was calibrated by fitting experimental breakthrough curves at various flow rates and titers



Performance objectives

- The model was used to numerically optimize the 1-4 column processes with respect to the following objectives:
- Productivity (Throughput), mAb produced per resin volume and time:




$$Prod = \frac{m_{mAb,eluate,cycle}}{n \cdot V_{col} \cdot t_{cycle}}$$

- Capacity Utilization (proportional to 1/resin costs): The Load, divided by static binding capacity. The higher the CU, the lower the resin costs.

$$CU = \frac{Load}{SBC} \sim \frac{1}{resin\ cost}$$

- Optimization variables: Flow rates, switch times

Optimization Results



2C-PCC outperforms other processes:

- Outperforms single column capture because of sequential loading principle, leading to increased capacity utilization (lower resin costs) and lower buffer demand
- Outperforms 3C-PCC and 4C-PCC in terms of productivity because its sequential loading phase is independent of recovery and regeneration steps running in parallel. Therefore the duration of sequential loading phase of 2C-PCC can be changed as larger titers are used

Conclusion

- Multicolumn process enable high capacity utilization and high throughput at the same time
- Multicolumn processes have 40-60% reduced resin costs, decreased buffer consumption and increased product concentration compared to batch chromatography

- 2C-PCC outperforms 3C-PCC and 4C-PCC in terms of productivity, while operating at similar capacity utilization and buffer consumption
- 2C-PCC requires least complex hardware of all multicolumn process, positive impact on equipment costs and risk of failure

Reference: Baur D et al. 2016. Comparison of batch and continuous multi-column protein A capture processes by optimal design. Biotechnology Journal. 11: 1860-7314